

Filament MVS Values

The following are MVS values calculated using the method described by [Muppet Labs](#). The value in the MVS column is the raw MVS, and should be reduced by about $\sim 0.5 \text{ mm}^3/\text{s}$.

Filament Brand	Filament Type	Nozzle Type	Nozzle Size	Temp (C)	Max Feedrate (mm/min)	Max Feedrate (mm/s)	MVS (mm^3/s)
Polaroid	PLA	Microswiss	0.6	200	485	8.083	19.44
Inland	PLA	Microswiss	0.6	225	605	10.083	24.25

Revision #3

Created 2021-02-14 15:45:07 UTC by Limb

Updated 2021-07-15 01:47:43 UTC by Limb